LC-MS User Guide

Release 1.0.0

Dominic Davis-Foster

May 17, 2020

Table of Contents

1	Befo	Before the Analysis									
	1.1	Sample Preparation									
	1.2	Pre-Run									
	1.3	Changing Eluent									
2	2 Mass Calibration / Check										
3	Perfe	orming Analysis									
	3.1	Starting a Run									
	3.2	During the Run									
	3.3	Shutdown									

CHAPTER 1

Before the Analysis

1.1 Sample Preparation

- Filter samples before injection with a Syringe and Filter.
 - Syringe filters can be found in the cupboard in the weighing room.
 - Clean syringes can be found in sealed packets in the cupboard under the HPLC.

1.2 Pre-Run

1. Check the indicator lights on the Autosampler, Pump, Column, VWD and MS are not Red.

```
Orange, Flashing Orange or Off are OK.
```

- 2. Check the valve near the ESI on the MS is open (handle in line with pipe) and is connected to L on the nearby block (shown below)
- 3. Check the levels of the bottles on the front of the MS. Tap the bottles to disturb the liquid if the level is not obvious.
- 4. Check the levels of the waste bottles in the following locations:
 - Beneath the bench below the Pump/Autosampler/Eluent stack
 - Beneath the bench below the MS and behind the Vacuum Pumps
 - Behind the MS
- 5. Check there is pressure in the Collision Gas Bottle behind the MS, and record the pressure
- 6. Check the 3 Exhaust Pipes are entering the Extractor for the Graphite Furnace
- 7. Return to the front of the instrument
- 8. Fetch some blue towel.



Fig. 1: Check that the pipe is connected to the port marked "L" to the left of this image, not to the port marked "H"

- 9. Open the Ionisation Chamber of the MS and spray 50% IPA around (but not into) the metal surfaces of the inlet orifice.
- 10. Spray the inside of the Ionisation Chamber door and wipe with the blue towel.
- 11. Close the Ionisation Chamber.
- 12. Open Mass Hunter Workstation Data Acquisition.
- 13. Check the pressures on the top right of the screen are within the following values:

Rough Vac	< 4.4×10 ⁰
Quad Vac	$1 \times 10^{5} - 4 \times 10^{5}$
TOC Vac	2×10 ⁷ - 4×10 ⁷

- 14. Check sufficient eluent is present in the eluent bottles being used for the run.
- 15. Check that the Quat. Pump is aware of the level of eluent in the bottles. To do this:
 - 1. Right click on the Quat. Pump panel of the *Instrument Status* window in *MassHunter*. <could use image here>
 - 2. Select *Bottle Fillings* from the menu.
 - 3. Under Actual Volume, ensure the value reflects the volume of eluent actually present in the bottle.
 - 4. Click OK.

Repeat the above steps when Changing Eluent (see below)

1.3 Changing Eluent

To change an eluent, perform the following steps:

- 1. Right click on the QTOF panel of the *Instrument Status* window in *MassHunter* and select $LC \rightarrow Waste$ from the menu.
- 2. In the Method Editor, under Quat. Pump, set the Flow Rate to 0.000 mL/min and press Enter

Bottle Fillings – 🗆 🔿												
Solvent Bottle												
Fillings												
	Actual Volume		Total Volume									
A:	0.28	liter	1.00 📫	liter								
B:	0.15	liter	1.00									
C:	0.44	liter	1.00	liter								
D:	0.38 🛟	liter	1.00 📫	liter								
Actions												
✓ Prevent analysis if level falls below 0.10 ÷												
✓ Turn pump off if running out of solvent												
			Ok	Cancel	Hel	P						

Fig. 2: The "Bottle Fillings" window



Fig. 3: The Instrument Status window

- 3. Swap the the eluent bottles. Try to swap similar eluents; for example, replace *Ammonium Formate* with *Ammonium Acetate*.
 - The Lids on the bottles do not need to be tight.
- 4. Open the valve behind the door on the front of the Pump.
- 5. Set the level of any solvents that weren't changed to 0% and enable the eluents that were changed.
- 6. Set the Flow Rate to 1.000 mL/min and press Enter
- 7. Check the pressure in the Column does not rise above 0 Bar. If it does, check the valve at the front of the Pump is open.
- 8. Set the Flow Rate to 5.000 mL/min and press Enter
- 9. After 5 minutes, set the Flow Rate to 0.000 mL/min and press Enter
- 10. Close the valve on the front of the Pump and close the door.

A method can be setup for flushing and equilibrating the column as part of the Worklist if desired

CHAPTER 2

Mass Calibration / Check

- 1. Select *Context* \rightarrow *Tune* at the top left of the *MassHunter* window.
- 2. Turn on the QTOF by clicking the On button in its section of the MassHunter window.
- 3. Close the valve near the ESI on the MS.
- 4. Set Calibrant Bottle to B at the bottom left of the Tune window.
- 5. Set LC Flow to to Waste at the bottom left of the Tune window.
- 6. In the Tune & Calibration tab (pictured below) set the instrument parameters required for the run, e.g.:
 - *Positive* if using +ESI
 - *TOF* <when would you use both?>
 - Mass Calibration / Check
 - 50-1500 m/z

From Agilent's manual:

```
Before you calibrate the instrument, you have to set the instrument state
to the proper instrument mode, mass range and fast polarity switching
mode. You set these values on the Instrument State tab.
When you change the mass range or enable/disable fast polarity switching
on the Instrument State tab, the pulser frequency is changed which results
in the DEI pulser warming up or cooling down. If the calibration is
performed too soon, the DEI may still be heating up or cooling down which
can result in drift. See the online Help for more information on the
Instrument State tab.
```

7. Click the button to start the calibration.

Note: The text of this button veries depending on the mode selected.

• For *Quadrupole*, the button reads *Start Check Quadrupole*

Tune File: TOFMassCalibration-1700mzRange.tun				Tune 8	Calibration Instrum	nent State					
Ion Polarity Positive C Negative 				Q-TDF: Low (1700 m/z) Extended Dynamic Range							
lon Source					 Positive Negative 	C G	Quadrupole TOF	 Mass Calibra Transmission 	tion / Check	 ⊙ 50-1700 m/z ⊂ 50-750 m/z 	
Dual AJS ESI	-				,		D			○ 50-250 m/z	
Gas Temp	225	243	°C			0	Both			Fragile Ions	
Drying Gas	13	12.9	1/min								
Nebulizer	20	20	psi						~		
VCap	4000 V	0.618	μΑ							Start TOF Mass	
Chamber		5.33	μΑ							Calibration	
Nozzle Voltage	2000 V										
Sheath Gas Tem	P 275	275	°C								
Sheath Gas Flow	12	12.0	l/min								
										Tune Report	
Calibrant Bottle	O None (0 A 🖲	В						~		
LC Flow to	● Waste (O MS		,							Apply

Fig. 1: The Tune & Calibration tab

- For TOF, the button reads Start TOF Mass Calibration
- For Both, the button reads Start Check Quadrupole & TOF Mass Calibration
- 8. Once complete, the calibration report will open automatically.
- 9. Under TOF Mass Calibration Data, for the largest mass check that the resolution is at least 20,000 and that the corrected residuals is below 5 ppm.
- 10. On Page 2 of the report, under *Detector* \rightarrow *MCP*, check that the value is below 900.
- 11. Back in *MassHunter*, return to *Acquisition* mode by selecting *Context* \rightarrow *Acquisition* at the top left of the window.
- 12. Open the valve near the ESI on the MS.

CHAPTER 3

Performing Analysis

3.1 Starting a Run

If Mass Calibration / Check was performed the LC flow will have been set to go to the MS automatically. Otherwise, perform the following steps to set the flow:

- 1. Right click on the QTOF panel of the Instrument Status window in MassHunter
- 2. Select $LC \rightarrow Waste$ from the menu.



Fig. 1: The Instrument Status window

To start the run, in the *Worklist*, select the samples to run by ticking the box to the left of the *Sample Name* field and click the button.

3.2 During the Run

1. Check that the Column pressure remains stable and does not rise above 200 bar. Analysis should stop automatically if the pressure rises above 600 bar.

i Worklist X											
10	7	📕 🚺 🕨 🔲 🛛 🖉	SulfaWorklist.wkl		•						
		Sample Name	Sample Position	Method	Data File	Sample Type	Level Name	Comment	Sample Group	Info. 📤	
1	v	Sulfa1	P1-A1	AJSESI_default.m	sulfas001.d	Sample					
2	V	Sulfa1	P1-A2	AJSESI_default.m	sulfas002.d	Sample					
3	v	Sulfa1	P1-A3	AJSESI_default.m	sulfas003.d	Sample				-	
4											
Worklist											
Worklist Method Editor Sample Run											

Fig. 2: The Worklist

2. Check that the three symbols in the Q-TOF panel of the *Instrument Status* window (shown below) remain green to indicate that the reference ions were found.



Fig. 3: The Instrument Status window

3.3 Shutdown

- 1. In the Instrument Status window, click the *Off* button to the lower right.
- 2. Close Data Acquisition.